

WATCO NuFit™ Tub Trim Kits Eliminate The Need To Remove Old Strainer Bodies!


How This New Patent Pending Design Works...



A. Remove the stopper from the old strainer body in the tub. Leave the old strainer body in the tub and place epoxy adhesive (furnished) on the flange.



B. Place the NuFit body into the old strainer body pressing the covering flange down onto the adhesive.




C. The NuFit body, with its bright new finish, is now securely in place ready for the new stopper. So quick and so easy!

NuFit Tub Trim Kits Are Available With Either A...



PRESFLO™ Stopper

Or



PUSH PULL™ Stopper

**YES! IT'S
ANOTHER INNOVATION
FROM WATCO!**



Always A Step Ahead

**JUST ANOTHER
EXAMPLE FROM THE
"ALWAYS A STEP AHEAD"
COMPANY!**

See back for more details